



Determination of Fibre

Fibre 6 system for fibre determination!

Fibre 6 apparatus are used in the determination of crude fiber in cereals, feed, food and other agricultural products.



The BEGER company established that the tedious analysis become fast, accurate and cost-effectively. We are glad to introduce the system for crude fibre determination – the FIBRE 6.

Fibre 6 apparatus are used in the determination of crude fiber in cereals, feed, food and other agricultural products. Due to its versatility and the ability of analysing up to 6 samples simultaneously, it becomes an indispensable tool in any laboratory that cannot be overlooked by the manufacturers of wide range of products required such an important parameter as the fibre content.

- ✓ easy in use
- ✓ compact size
- ✓ affordability

Specifications:

Number of samples	6
Set of filter bags	30
Rated voltage	230 V / 50 Hz
Power consumption	500 W
Dimensions (W x H x D)	155 x 580 x 235 mm
Weight	4 kg
Heat control	Power gain
Cat. №	314 106 00 00



**Stegne 27, Industrijska cona Stegne,
1000 Ljubljana, Slovenia**

www.beger.si

e-mail: info@beger.si

skype: [beger_d.o.o](https://www.skype.com/people/beger_d.o.o)

Find us on:

