

<p><b>Vertical Pressure Steam Sterilizer</b></p>	<p>Valid volume: 120L(φ 480*660mm)            150L(φ 510*740mm)            Working pressure: 0.22Mpa            Working temperature:134°C            Max. working pressure:0.23Mpa            Heat average: ≤ ±1°C            Timer: 0-60min            Adjustment of temperature:105-134°C            Power: 6kw/AC220V.50Hz            Transport dimension:730*710*1250(mm)            760*760*1280mm            G/N.W:120kg/100kg 130kg/105kg</p>		<p>With timing, temperature control, below exhaust, water cut</p>
<p><b>Vertical Pressure Steam Sterilizer</b></p>	<p>Valid volume:35L(φ 318*450mm)            50L(φ 340*550mm) 75L(φ 400*600mm)            100L(φ 440*650)            Working pressure:0.22Mpa            Working temperature:134°C            Max. working pressure:0.23Mpa            Heat average: ≤ ±1°C            Timer: 0-60min            Adjustment of temperature:105-134°C            Power:2.5kw/AC220V.50Hz 3kw/AC220V.50Hz            4.5kw/AC220V.50Hz 4.5kw/AC220.50hz            Transport dimension:590*550*880(mm)            590*590*1120(mm)            650*620*1120(mm)670*650*1200(mm)            G/N.W:55kg/45kg 70kg/50kg 90kg/69kg            100kg/76kg</p>		<p>With timing, temperature adjuster, water cut</p>
<p><b>Vertical Pressure Steam Sterilizer (hand round automatic)</b></p>	<p>Valid volume:35L(φ 318*450mm)            50L(φ 340*550mm) 75L(φ 400*600mm)            100l(φ 440*650)            Working pressure:0.22Mpa 0.22mpa 0.22mpa            0.14mpa            Working temperature:134°C 134°C 134°C 126°C            Max. working pressure:0.23Mpa0.23Mpa0.23Mpa0.165Mpa            Heat average: ≤ ±1°C            Timer: 0-60min            Adjustment of temperature:105-134°C 105-134°C105-134°C105-126°C            Power:2.5kw/AC220V.50Hz 3kw/AC220V.50Hz            4.5kw/AC220V.50Hz 4.5kw/AC220.50hz            Transport dimension:510*480*1100(mm)            540*510*1150(mm) 650*630*1300(mm)            G/N.W:62kg/50kg 65kg/55kg 105kg/82kg            105kg/95kg</p>		<p>Hand round automatic, with timing, temperature control, water cut</p>




<p><b>Vertical pressure steam sterilizer (fully automatic autoclave)</b></p>	<p>Valid volume:35L(φ 318*450mm) 50L(φ 340*550mm) 75L(φ 400*600mm) 100L(440*650mm)</p> <p>Working pressure:0.22Mpa0.22mpa0.22mpa0.14mpa</p> <p>Working temperature:134°C134°C134°C126°C</p> <p>Max. working pressure:0.23Mpa0.23Mpa0.23Mpa0.165mpa</p> <p>Heat average: ≤ ±1°C</p> <p>Timer: 0-60min</p> <p>Adjustment of temperature:105-134°C105-134°C105-134°C 105-126°C</p> <p>Power:2.5kw/AC220V.50Hz 3kw/AC220V.50Hz 4.5kw/AC220V.50Hz 4.5kw/AC220.50hz</p> <p>Transport dimension:510*480*1100(mm) 540*510*1150(mm) 650*630*1300(mm) 630*650*1400mm</p> <p>G/N.W:62kg/50kg 65kg/55kg 105kg/82kg 105kg/95kg</p>		<p>Hand round automatic, digital display, disconnect air conditioning, water cut, automatically disconnect the steam</p>
<p><b>Vertical pressure steam sterilizer (hand round automatic)</b></p>	<p>Valid volume:35L(φ 318*450mm) 50L(φ 340*550mm) 75L(φ 400*600mm) 100L(φ 440*650mm)</p> <p>Working pressure:0.22Mpa0.22mpa0.22mpa0.14mpa</p> <p>Working temperature:134°C134°C134°C126°C</p> <p>Max. working pressure:0.23Mpa0.23Mpa0.23Mpa0.165Mpa</p> <p>Heat average: ≤ ±1°C</p> <p>Timer: 0-60min</p> <p>Adjustment of temperature:105-134°C105-134°C105-134°C105-126°C</p> <p>Power:2.5kw/AC220V.50Hz 3kw/AC220V.50Hz 4.5kw/AC220V.50Hz 4.5kw/AC220.50hz</p> <p>Transport dimension:500*650*1100(mm) 500*670*1150(mm) 550*750*1260(mm)620*820*1400mm</p> <p>G/N.W:62kg/50kg 80kg/55kg 105kg/82kg 134kg/100kg</p>		<p>With timing, temperature control, water cut</p>

<p><b>Vertical pressure steam sterilizer (fully automatic microcomputer)</b></p>	<p>Valid volume:35L(φ 318*450mm) 50L(φ 340*550mm) 75L(φ 400*600mm) 100(φ 440*650mm)</p> <p>Working pressure:0.22Mpa Working temperature:134°C Max. working pressure:0.23Mpa Heat average: ≤ ±1°C Timer: 0-99minOR 0-99hour59min Adjustment of temperature:105-134°C Power:2.5kw/AC220V.50Hz 3kw/AC220V.50Hz 4.5kw/AC220V.50Hz 4.5kw/AC220.50hz Transport dimension:550*550*880(mm) 590*590*1120(mm) 650*620*1120(mm)670*650*1200mm</p> <p>G/N.W:55kg/45kg 70kg/50kg 90kg/69kg 100kg/76kg</p>		<p>Automatically digital display, disconnect air conditioning, water cut, automatically disconnect the steam</p>
<p><b>Vertical pressure steam sterilizer (fully automatic microcomputer)</b></p>	<p>Valid volume:35L(φ 318*450mm) 50L(φ 340*550mm) 75L(φ 400*600mm ) 100l(φ 440*650)</p> <p>Working pressure:0.22Mpa0.22mpa0.22mpa0.14mpa</p> <p>Working temperature:134°C 134°C 134°C 126°C Max. working pressure:0.23Mpa 0.23Mpa 0.23Mpa 0.165Mpa Heat average: ≤ ±1°C Timer: 0-99min or0~99hour 59min Adjustment of temperature:0-134°C 0-134°C0-134°C0-126°C</p> <p>Power:2.5kw/AC220V.50Hz 3kw/AC220V.50Hz 4.5kw/AC220V.50Hz 4.5kw/AC220.50hz</p> <p>Transport dimension:500*650*1100(mm) 500*690*1150(mm) 550*750*1260(mm) 620*820*1400mm</p> <p>G/N.W:62kg/50kg 65kg/55kg 105kg/82kg 115kg100kg</p>		<p>Digital display, disconnect air conditioning, water cut,automatically disconnect the steam</p>
<p><b>Vertical pressure steam sterilizer</b></p>	<p>Valid volume:35L(φ 318*450mm) Working temperature:126°C Max. working pressure:0~0.16Mpa Timer: 0-60min Power:2.5kw/AC220V.50Hz Transport dimension:550*550*880(mm) G/N.W:55kg/45kg</p>		<p>Outer covering carbon steel power coating, with timing</p>

<p><b>Vertical pressure steam sterilizer</b></p>	<p>Valid volume:50L(φ 340*550mm)                  Working temperature:126°C                  Max. working pressure:0.14Mpa                  Timer: 0-60min                  Power:4.5kw/AC220V.50Hz                  Transport dimension:590*550*1100(mm)                  G/N.W:65kg/55kg</p>		<p>Outer covering carbon steel power coating, with timing</p>
<p><b>Table top steam sterilizer</b></p>	<p>Valid volume:20L(φ 270*340mm)                  Max. working pressure:0.22Mpa                  Max.working temperature:134°C                  Adjustment of temperature:105-134°C                  Heat average: ≤ ±1°C                  Power:1.5kw/AC220V.50Hz                  Sterilizing net:φ 235*255mm                  Transport dimension:420*420*630(mm)                  G/N.W:29/27kg</p>		<p>Max. working temperature:121</p>
<p><b>Table top steam sterilizer</b></p>	<p>Valid volume:20L(φ 250*420mm)                  25L(φ 250*520mm)                  Max. working pressure:0.22Mpa                  Max.working temperature:134°C                  Adjustment of temperature:105-134°C                  Timer: 0-60min                  Heat average: ≤ ±1°C                  Power:1.5kw/AC220V.50Hz (AC110V 60Hz)                  Sterilizing plate:340*200*30mm(3piece)                  400*200*30mm(3 piece)                  Transport dimension:700*580*500(mm)                  800*580*500mm                  G/N.W : 43/38kg 50/44kg</p>		<p>N class</p>
<p><b>Table top steam sterilizer (fully automatic microcomputer type)</b></p>	<p>Valid volume:20L(φ 250*420mm)                  25L(φ 250*520mm)                  Max. working pressure:0.22Mpa                  Max.working temperature:134°C                  Adjustment of temperature:105-134°C                  Timer: 0-99min                  Heat average: ≤ ±1°C                  Power:1.5kw/AC220V.50Hz                  Sterilizing plate:340*200*30mm(3piece)                  400*200*30mm(3 piece)                  400*200*30mm(3 pieces)                  G/N.W : 43/38kg 50/45kg</p>		<p>Computer control, N class</p>

<p><b>Table top steam sterilizer</b></p>	<p>Valid volume: 20L(φ 250*420mm) 25L(φ 250*520mm) Max. working pressure:0.22Mpa Max.working temperature:134°C Adjustment of temperature:105-134°C Timer: 0-60min Heat average: ≤ ±1°C Power:1.5kw/AC220V.50Hz (AC110V 50Hz) Sterilizing plate:340*200*30mm(3piece) 400*200*30mm(3 piece)  Transportdimension:700*580*500(mm) 800*580*500mm  G/N.W:44/40kg 50/45kg</p>		<p>N class</p>
<p><b>Table top steam sterilizer</b></p>	<p>Valid volume:35L(φ 300*500mm) 50L(φ 340*550mm) Max. working pressure:0.22Mpa Max.working temperature:134°C Adjustment of temperature:105-134°C Timer: 0-60min Heat average: ≤ ±1°C Power:2.7kw/AC220V.50Hz (AC110V 60Hz) Sterilizing plate:400*200*30mm(4piece) 500*250*30mm(4 piece) G/N.W : 95/80kg 93/110kg</p>		<p>N class</p>
<p><b>Table top steam sterilizer (fully automatic microcomputer type)</b></p>	<p>Valid volume:35L(φ 300*500mm) 50L(φ 340*550mm) Max. working pressure:0.22Mpa Max.working temperature:134°C Temp.range: 105 ~ 134°C Timer: 0 ~ 99min Power:2.7kw/AC220V.50Hz Sterilizing plate:400*200*30mm(4piece) 500*250*30mm(4 piece) Dimension: 73*55*51 ( cm ) 77*60*54cm G/N.W : 94/80kg 110kg/92kg</p>		<p>Computer control, N class</p>

<p><b>Table top steam sterilizer</b></p>	<p>Valid volume:20L(φ 270*3400mm)  Max. working pressure:0.22Mpa  Max.working temperature:134°C  Adjustment of temperature:105-134°C  Timer: 0-60min  Heat average: ≤ ±1°C  Power:1.5kw/AC220V.50Hz  Sterilizing net:φ 235*255mm  Transport dimension:420*420*630(mm)  G/N.W : 29/27kg</p>		<p>N class</p>
<p><b>Table top steam sterilizer</b></p>	<p>Valid volume:25L(φ 270*410mm)  Max. working pressure:0.22Mpa  Max.working temperature:134°C  Adjustment of temperature:105-134°C  Timer: 0-60min  Heat average: ≤ ±1°C  Power:1.5kw/AC220V.50Hz  Sterilizing plate:340*200*30mm(3piece)  Transport dimension:620*410*570(mm)  G/N.W : 19/18kg</p>		<p>N class</p>
<p><b>Portable pressure steam sterilizer (Electric or LPG heated)</b></p>	<p>Sterilizer volume:12L(φ 240*200mm)  Working pressure:0.14-0.16Mpa  Working temperature:126°C  Max. safety pressure:0.165Mpa  Power:AC220V.50Hz/2KW  Dimension:360*360*420(mm)  N/G.W:11/12KG</p>		<p>Electric or LPG</p>
<p><b>Portable pressure steam sterilizer</b></p>	<p>Sterilizer volume:18L(φ 280*260mm)  24L(φ 280*390mm)  Working pressure:0.14-0.16Mpa  Working temperature:126°C  Max. safety pressure:0.165Mpa  Heat average: ≤ ±1°C  Power:AC220V.50Hz/2KW  Dimension:410*410*520(mm) 410*410*630(mm)  N/G.W16/18KG 17/19KG</p>		

<p><b>Portable pressure steam sterilizer</b></p>	<p>Sterilizer volume:18L(φ 280*260mm) 24L(φ 290*380)                  Working pressure:0.14~0.16Mpa                  Working temperature:126°C                  Max. safety pressure:0.165Mpa                  Power:AC220V.50Hz/2KW                  Dimension:410*410*430(mm) 410*410*630mm                  N/G.W:14/16KG 15/17KG</p>		<p>Electric or LPG heated</p>
<p><b>Portable pressure steam sterilizer</b></p>	<p>Sterilizer volume:18L(φ 280*265mm)                  Working pressure:0.14-0.16Mpa                  Working temperature:126°C                  Max. safety pressure:0.165Mpa                  Power:AC220V.50Hz/2KW                  Dimension:410*410*540(mm)                  N/G.W:11/13KG</p>		<p>Electric or LPG heated</p>
<p><b>Portable pressure steam sterilizer</b></p>	<p>Sterilizer volume:18L(φ 280*260mm)                  24L(φ 280*390mm)                  Working pressure:0.14~0.16Mpa                  Working temperature:126°C                  Max. safety pressure:0.165Mpa                  Heat average: ≤ ±1°C                  Power:AC220V.50Hz/2KW                  Dimension:410*410*630(mm) 410*410*650(mm)                  N/G.W:12/13KG 15/14KG</p>		<p>With timing, temperature adjuster, volume:18L, 24L</p>
<p><b>Portable steam sterilizer (Automatic microcomputer type)</b></p>	<p>Sterilizer volume:18L(φ 280*260mm)                  Working pressure:0.14-0.16Mpa                  Working temperature:126°C                  Temperature adjust range:105-126°C                  Timer range:0-99min                  Max. safety pressure:0.165Mpa                  Heat average: ≤ ±1°C                  Power:AC220V.50Hz/2KW                  Dimension:410*410*630(mm)                  N/G.W:12/14KG</p>		<p>microcomputer control, volume:18L</p>

<p><b>Portable steam sterilizer</b></p>	<p>Sterilizer volume:18L(φ 280*260mm)                  24L(φ 280*390mm)                  Working pressure:0.14~0.16Mpa                  Working temperature:126°C                  Temperature adjust range:105-126°C                  Timer range:0-60min                  Max. safety pressure:0.165Mpa                  Heat average: ≤ ±1°C                  Power:AC220V.50Hz/2KW                  Dimension:410*410*630(mm) 410*410*750(mm)                  N/G.W:12/13KG 14/15KG</p>		<p>Water cut,                  volume:18L, 24L</p>
<p><b>Wall mounted water distillation electric</b></p>	<p>Power : AC220V~240V/50Hz~60Hz                  Output : 3.5~4.5L/h                  Consumption : 1.5Kw*2                  Net weight : 7kg/set                  Gross weight : 8kg/set 17kg/2set                  Dimension : 740*400*140mm/set                  760*420*310mm/2set                  Capacity of pan : 5.5liter PH:5~6</p>		<p>Wall mounted</p>
<p><b>Water distilling apparatus tower type electrically heated</b></p>	<p>Output :5liters/hour                  10liters/hour                  20liters/hour                  Power : AC220v/4.5Kw                  AC380v/7.5Kw                  AC380v/13.5Kw                  Package : 540*270*490mm                  610*310*570mm                  740*400*680mm                  Weight : 9kg 11kg 17kg</p>		<p>Tower type</p>